



Intramesh SOFT LIFT

Suburethral support tape

The **easy** way




Photos and texts non-contractual. Specifications can be modified without notice.
Cousin Biotech S.A.S. with a capital of 340 656 € - 398 460 261 RCS Lille - N°TVA FR 34 398 460 261

Cousin Biotech
8, rue de l'Abbé Bonpain
59117 Wervicq Sud, France
Tel.: +33 (0) 3 20 14 40 00
Fax: +33 (0) 3 20 14 41 22
mail: uro-gyneco@cousin-biotech.com
www.cousin-biotech.com



Distributed by:




Intramesh SOFT LIFT is a class IIIb medical device manufactured by Cousin Biotech s.a.s.
Please read instructions for use carefully.
The management system of COUSIN BIOTECH s.a.s has been certified as meeting the requirements of ISO 13485.
Reference: PLGSLG03 - 01/07/13

Urology and gynaecology surgery


Intramesh SOFT LIFT

Solution for treatment of stress urinary incontinence in women


The easy way



Polypropylene
Monofilament
Macroporous



Winged guide
for TOT In-Out




Single use ancillaries

Material

- Knitted polypropylene, monofilament
- Knitted in its final width
- Macroporous
- Shape memory
- Resistant

Description

Suburethral support tape for treatment of stress urinary incontinence in women. Its unique knit ensures an ideal positioning of the sling under the urethra



Advantages

New splice-gripping system allowing a smooth connection and disconnection of the implant from the needle in a single move

Safety

- Well tolerated material
- No particles released
- Optimal recolonization
- No « rope effect »
- Minimal shrinkage

Patient's comfort

- Day surgery
- Less implanted material
- Edges of the tape are not rough
- Less traumatic procedure


Stock management

- One kit = 4 surgical approaches
- Complete range of ancillaries allowing TVT and TOT approaches
- Winged guide included for TOT In-Out
- Single use kit

References

Reference	Description and weight
GCB INSOFT U	Single use kit SOFT LIFT for stress urinary incontinence, 60 g/m ²
GCB INLIFT U	Single use kit LIFT for stress urinary incontinence, 100 g/m ²

Single use kits include

1 implant	2 transobturator ancillaries	1 retropubic ancillary	1 guide
 (1 cm x 45 cm)	